


# geneMAG-RNA/DNA 96

compatible for KingFisher™ 96 and BioSprint™ 96

Magnetic RNA/DNA purification kit  
compatible for KingFisher 96™ or  
BioSprint 96™ workstation.




**chemicell**

NEW TOOLS IN BIOSCIENCES

# Technology

The **geneMAG-RNA/DNA 96** kit is a novel, simple and highly efficient tool for the isolation of total RNA/DNA with magnetic silica beads.

Simple washing steps with two different buffers remove salts, metabolites and macromolecular cellular components. Pure RNA/DNA is finally eluted under low ionic strength conditions with RNase-free water.



The magnetic bead processing of KingFisher™ 96 / BioSprint™ 96 workstation

The KingFisher™ 96 workstation is a trademark of Thermo Fisher Scientific.

The BioSprint™ 96 workstation is a trademark of Qiagen.

# Products

Kits	Contents	Number of isolations	Price Euro/US\$
<b>geneMAG-RNA/DNA 96</b> (Cat. No.: KF3401-96)	<ul style="list-style-type: none"><li>• 100 ml Lysis &amp; Binding Buffer</li><li>• 200 ml Wash Buffer I</li><li>• 10 ml SiMAG/KF-DNA Beads</li></ul>	1 x 96 preps per 10 <sup>9</sup> bacteria	185 / 240
<b>geneMAG-RNA/DNA 480</b> (Cat. No.: KF3401-480)	<ul style="list-style-type: none"><li>• 500 ml Lysis &amp; Binding Buffer</li><li>• 1000 ml Wash Buffer I</li><li>• 50 ml SiMAG/KF-DNA Beads</li></ul>	5 x 96 preps per 10 <sup>9</sup> bacteria	725 / 942

## Reagents and Equipment to be Supplied by the User

- **Wash Buffer II:** 70% Ethanol or 70% Isopropanol.
- **Elution Buffer:** Nuclease-free water or DEPC-Water for elution of RNA/DNA from the beads.
- **KingFisher™ 96 / BioSprint™ 96 workstation**
- **Deep well 96-well plates (2,2 ml) squared well**
- **KingFisher™ 96 plate (0,3 ml)**
- **Magnet Head for deep well 96-well plates**
- **DNase Treatment:** DNA-free™ Kit, DNase Treatment and Removal Reagents Part Number AM1906 (Applied Biosystems).

## Safety Note

**Lysis/Binding-Buffer** and **Wash Buffer I** contain chaotropic salts, which are irritant. Take appropriate laboratory safety measures and wear gloves when handling. **Avoid skin and eye contact**

## Protocol for KingFisher™ 96 or BioSprint 96™

This protocol describes the isolation of genomic-DNA from  $10^9$  bacteria cells per each well of 96-well plate with the geneMAG-RNA/DNA 96 kit using KingFisher™ 96 or BioSprint™ 96 workstation.

**Preparation of the deep well 96-well plates (2,2 ml)**

### **First 96-well plate:**

1. Add 100 µl **Bacteria Suspension** in each well

**Bacteria Suspension:** Add 1.5 ml of overnight cultured cells (approximately  $10^8$  cells) into a 1.5 ml microcentrifuge tube. Centrifuge for 2 minutes at 11,000 x g to pellet the cells. Discard the supernatant.

Resuspend the bacteria pellet in 100 µl **Lysis & Binding Buffer**

2. Add 1000 µl **Lysis & Binding Buffer** and 100 µl **SiMAG/KF-DNA**.

### **Second 96-well plate:**

1. Add 1000 µl **Wash Buffer I** in each well

### **Third 96-well plate:**

1. Add 1000 µl **Wash Buffer I** in each well

### **Fourth 96-well plate:**

1. Add 1000 µl **Wash Buffer II** (70% 2-propanol) in each well

### **Fifth 96-well plate:**

1. Add 1000 µl **Wash Buffer II** (70% 2-propanol) in each well

### **Sixth 96-well plate: Use 96-well plate with max. volume of 0,3 ml**

1. Add 100 µl **Elution Buffer** (dH<sub>2</sub>O) in each well


### **Seventh 96-well plate: Parking station!!!!**

## Protocol for KingFisher™ 96 or BioSprint 96™

**Settings of the processing times for Lysis- & Binding-, Wash- and Elution steps :**

**Start the KingFisher Software and set the following parameters:**

1. Lysis & Binding process (first deep-well plate):  
6 minutes with low stirring
2. Wash process with Wash Buffer I (second deep-well plate):  
1 minutes with medium stirring
3. Wash process with Wash Buffer I (third deep-well plate):  
1 minutes with medium stirring
4. Wash process with Wash Buffer II (fourth deep-well plate):  
1 minutes with medium stirring
5. Wash process with Wash Buffer II (fifth deep-well plate):  
1 minutes with medium stirring
6. Elution process with Elution Buffer (dH<sub>2</sub>O):  
Heat time: 10 minutes with high stirring  
Temperature: 80°C




Total RNA/DNA was isolated from 1.5ml E. coli LB culture using geneMAG-RNA/DNA 96.  
(Data kindly provided by Cengiz Öztürk, Charité, University Hospital of Humboldt-University to Berlin, Germany)

## Protocol

### Optional: DNase Treatment

**Use Ambion® DNA-free™ DNase Treatment and Removal Reagents and follow the instructions of the manufacturer.**

10. Add 5 µL (2 Units/ µl) rDNase I to the eluted RNA/DNA and mix gently.
11. Incubate at 37°C for 30 min.
12. Add 10 µl resuspended DNase Inactivation Reagent and mix well.
13. Incubate 2 minutes at room temperature, mix occasionally.
14. Centrifuge at 11,000 x g for 2 minutes, carefully transfer the RNA containing supernatant into a fresh tube.



Total RNA was isolated from 1 ml E. coli LB culture using geneMAG-RNA/DNA 96 with subsequent DNase treatment (Ambion).

(Data kindly provided by Cengiz Öztürk, Charité, University Hospital of Humboldt-University to Berlin, Germany)

## Contact

### chemicell GmbH

Eresburgstrasse 22-23  
12103 Berlin  
Germany



[info@chemicell.com](mailto:info@chemicell.com)

Tel.: +49-30-2141481  
Fax.: +49-30-21913737  
e-mail: [info@chemicell.com](mailto:info@chemicell.com)  
Internet: [www.chemicell.com](http://www.chemicell.com)