
PEG Product Information

PEG-di-Cibacron Blue

| | |
|---------------------|---|
| Product: | α,ω-di-Cibacron Blue polyethylene glycol |
| Article Number: | PDC20-1 (1 g); PDC20-5 (5 g) |
| Description: | Aqueous solution of α,ω -di-Cibacron Blue polyethylene glycol |
| Weight of Volume: | - |
| Lot: | - |
| Production Date: | - |
| PEG _{MW} : | 20000 Da |
| Functional Group: | Cibacron Blue |
| Storage Buffer: | ddH ₂ O |
| Autoclaved: | Yes |
| Storage: | 4 – 8 °C. Protect from light; Do not freeze! |
| Expiry date: | Two years after production date |

