
PEG Product Information

PEG-mono-Cibacron Blue

Product:	α-Methoxy-ω-Cibacron Blue polyethylene glycol
Article Number:	PDC5-1 (1 g); PDC5-5 (5 g)
Description:	Aqueous solution of α -Methoxy- ω -Cibacron Blue polyethylene glycol
Weight of Volume:	-
Lot:	-
Production Date:	-
PEG _{MW} :	5000 Da
Functional Group:	Cibacron Blue
Storage Buffer:	ddH ₂ O
Autoclaved:	Yes
Storage:	4 – 8 °C. Protect from light; Do not freeze!
Expiry date:	Two years after production date

