

---

## PEG Product Information

### PEG-di-Cibacron Blue

---

|                     |   |
|---------------------|---|
| Product:            | <b><math>\alpha,\omega</math>-di-Cibacron Blue polyethylene glycol</b>    |
| Article Number:     | PDC6-1 (1 g); PDC6-5 (5 g)  |
| Description:        | Aqueous solution of $\alpha,\omega$ -di-Cibacron Blue polyethylene glycol |
| Weight of Volume:   | -   |
| Lot:                | -   |
| Production Date:    | -   |
| PEG <sub>MW</sub> : | 6000 Da   |
| Functional Group:   | Cibacron Blue   |
| Storage Buffer:     | ddH <sub>2</sub> O  |
| Autoclaved:         | Yes   |
| Storage:            | 4 – 8 °C. <b>Protect from light; Do not freeze!</b>                       |
| Expiry date:        | Two years after production date   |

---

