
Product Information – screenCORE-Goat anti-mouse IgG

Product:	screenCORE-Goat anti-mouse IgG (whole molecule) (Fluorescent Beads)				
Article Number:	6202-1 (1 ml); 6202-5 (5 ml)				
Description:	Aqueous dispersion of fluorescent silica particles				
Application:	For binding of secondary antibody				
Lot Number:					
Production Date:					
Weight of Volume:	10 mg/ml				
Matrix:	Silica, high-porous				
Size (hydrodynamic diameter):	1.0 µm				
Number of Particles:	~ 1.8 x 10 ¹² /g				
Surface Area:	~ 100 m ² /g				
Density:	~ 2.25 g/cm ³				
Functional Group:	Goat anti-mouse IgG (whole molecule)				
screenCORE/ Fluorescence Color:	B blue	GQ green	GF green	GA green	G green
Excitation:	400 nm	423 nm	436 nm	460 nm	502 nm
Emission:	420 nm	503 nm	520 nm	550 nm	525 nm
screenCORE/ Fluorescence Color:	OE orange	O orange	OP orange	RR red	R red
Excitation:	510 nm	526 nm	536 nm	540 nm	633 nm
Emission:	595 nm	555 nm	617 nm	625 nm	672 nm
Autoclaved:	No				
Storage Buffer / Solution:	PBS, 0.05 % sodium azide				
Storage:	At 4 – 8 °C. Do not freeze! PROTECT FROM LIGHT!				
Expiry Date:	Six months after production date.				
Note:	For complete resuspension vortex thoroughly!				



NOTE: The fluorescence of the screenMAG particles is only detectable on the same side where the excitation takes place.

Please note that there is a difference in fluorescence observation between dissolved fluorescence molecules and solid fluorescence particles. Fluorescence spectrophotometer with a fluorescence detection unit with an angle of 90° to the excitation source will detect no or only weak fluorescence signals.