
Product Information – fluidMAG-PS

Product:	fluidMAG-PS		
Article Number:	4111-1 (1 ml) ; 4111-5 (5 ml)		
Description:	Aqueous dispersion of magnetic nanoparticles		
Application:	Strong cation exchanger		
Weight of Volume:	25 mg/ml		
Production Date:			
Lot:			
Core:	Magnetite		
Matrix:	Poly(sodium 4-styrene-sulfonate)		
Size approx.:	50 nm	100 nm	200 nm
Number of Particles:	$\sim 1.3 \times 10^{16}/g$	$\sim 1.8 \times 10^{15}/g$	$\sim 2.2 \times 10^{14}/g$
Density:	$\sim 1.25 \text{ g/cm}^3$		
Type of Magnetization:	Superparamagnetic		
Functional Group:	Sodium Sulfonate, $-\text{SO}_3^- \text{Na}^+$		
Autoclaved:	Yes		
Storage Buffer:	ddH ₂ O		
Storage:	4 – 8 °C. Do not freeze!		
Expiry date:	Two years after production date		

