

---

## PEG Product Information

### PEG-di-Cibacron Blue

---

Product:	<b><math>\alpha,\omega</math>-di-Cibacron Blue polyethylene glycol</b>
Article Number:	PDC8-1 (1 g); PDC8-5 (5 g)
Description:	Aqueous solution of $\alpha,\omega$ -di-Cibacron Blue polyethylene glycol
Weight of Volume:	-
Lot:	-
Production Date:	-
PEG <sub>MW</sub> :	8000 Da
Functional Group:	Cibacron Blue
Storage Buffer:	ddH <sub>2</sub> O
Autoclaved:	Yes
Storage:	4 – 8 °C. <b>Protect from light; Do not freeze!</b>
Expiry date:	Two years after production date

---

