
PEG Product Information

PEG-di-Phosphate

Product:	α,ω-Di-phosphate polyethylene glycol
Article Number:	PDP6-1 (1 g); PDP6-5 (5 g)
Description:	Aqueous solution of α,ω -Di-phosphate polyethylene glycol
Weight of Volume:	-
Lot:	-
Production Date:	-
PEG _{MW} :	6000 Da
Functional Group:	Phosphate
Storage Buffer:	ddH ₂ O
Autoclaved:	Yes
Storage:	4 – 8 °C.
Expiry date:	Two years after production date

