

---

## PEG Product Information

### PEG-di-Succinic acid

---

Product:	<b><math>\alpha,\omega</math>-Di-succinic acid polyethylene glycol</b>
Article Number:	PDS20 (1 g); PDS20 (5 g)
Description:	Aqueous solution of $\alpha,\omega$ -Di-succinic acid poly-ethylene glycol
Weight of Volume:	-
Lot:	-
Production Date:	-
PEG <sub>MW</sub> :	20000 Da
Functional Group:	Carboxylate
Storage Buffer:	ddH <sub>2</sub> O
Autoclaved:	Yes
Storage:	4 – 8 °C.
Expiry date:	Two years after production date

---

