
PEG Product Information

PEG-di-Succinic acid

Product:	α,ω-Di-succinic acid polyethylene glycol
Article Number:	PDS6 (1 g); PDS6 (5 g)
Description:	Aqueous solution of α,ω -Di-succinic acid poly-ethylene glycol
Weight of Volume:	mg/ml
Lot:	
Production Date:	
PEG _{MW} :	6000 Da
Functional Group:	Carboxylate
Storage Buffer:	ddH ₂ O
Autoclaved:	Yes
Storage:	4 – 8 °C.
Expiry date:	Two years after production date

