
PEG Product Information

PEG-Di-ester benzene-dicarboxylate

Product:	α,ω-Di-1-ester benzene-2,4-dicarboxylate-polyethylene glycol
Article Number:	PDT8-1 (1 g); PDT8-5 (5 g)
Description:	Aqueous solution of α,ω -Di-1-ester benzene-2,4-dicarboxylate-polyethylene glycol
Weight of Volume:	-
Lot:	-
Production Date:	-
PEG _{MW} :	8000 Da
Functional Group:	Carboxylate
Storage Buffer:	ddH ₂ O
Autoclaved:	Yes
Storage:	4 – 8 °C. Protect from light; Do not freeze!
Expiry date:	Two years after production date

