
PEG Product Information

PEG-Mono-ester benzene-dicarboxylate

Product:	α-methoxy-ω-1-ester benzene-2,4-dicarboxylate polyethylene glycol
Article Number:	PMT5-1(1g); PMT5-5 (5 g)
Description:	Aqueous solution of α -methoxy- ω -1-ester benzene-2,4-dicarboxylate polyethylene glycol
Weight of Volume:	mg/ml
Lot:	
Production Date:	
PEG _{MW} :	5000 Da
Functional Group:	Carboxylate
Storage Buffer:	ddH ₂ O
Autoclaved:	Yes
Storage:	4 – 8 °C. Protect from light; Do not freeze!
Expiry date:	Two years after production date
